

BLACKJACK SPEC SHEET



- IMPROVE HCAHPS SCORES
- INCREASE NURSE CALL CONNECTIVITY
- DECREASE DAMAGE TO NURSE CALL SYSTEMS AND BEDS CAUSED BY NURSE CALL CABLES
- SAVE MONEY ON REPLACEMENT CABLES AND TECHNICIAN TIME
- UNIVERSALLY APPLICABLE TO ALL BED TYPES AND NURSE CALL SYSTEMS
- INDIVIDUAL SPRING LOADED, GOLD PLATED AND MAGNETICALLY BACK PINS
- CORROSION RESISTANT



Break	Bed Type	Product Photo	Part Number	Description
Wall Break	Hill-Rom Linet		XFP	Short 1/4" phono jack to magnetic end cable for wall + Long magnetic end to 37-pin conventional female cable for bed
	ArjoHuntleigh Stryker Unmano		XMP	Short 1/4" phono jack to magnetic end cable for wall + Long magnetic end to 37-pin conventional male cable for bed
	Hill-Rom Linet		XF6	Magnetic 37-pin male brick for wall + Long 37-pin magnetic female to 37-pin conventional female cable for bed
	ArjoHuntleigh Stryker Unmano		XM6	Magnetic 37-pin male brick for wall + Long 37-pin magnetic female to 37-pin conventional male cable for bed
	Hill-Rom Linet		XFR	Magnetic 37-pin male brick for wall + Long 37-pin magnetic right-angle female to 37-pin conventional female cable for bed
	ArjoHuntleigh Stryker Unmano		XMR	Magnetic 37-pin male brick for wall + Long 37-pin magnetic right-angle female to 37-pin conventional male cable for bed
	Hill-Rom Linet		XF6CFR2KIT	Short 37-pin conventional right-angle male to 37-pin magnetic male cable for wall + Long 37-pin magnetic female to 37-pin conventional female cable for bed
	ArjoHuntleigh Stryker Unmano		XM6CFR2KIT	Short 37-pin conventional right-angle male to 37-pin magnetic male cable for wall + Long 37-pin magnetic female to 37-pin conventional male cable for bed
Bed Break	Hill-Rom Linet		XF2CFRKIT	Long 37-pin conventional right-angle male to 37-pin magnetic male cable for wall + Short 37-pin magnetic female to 37-pin conventional female cable for bed
	ArjoHuntleigh Stryker Unmano		XM2CFRKIT	Long 37-pin conventional right-angle male to 37-pin magnetic male cable for wall + Short 37-pin magnetic female to 37-pin conventional male cable for bed